

NEW

InterActive™

Conical Connection with Added Simplicity



Conical Connection
compatible with NobelActive™*
& NobelReplace™ CC*

Square top detaches with
impression for metal to metal
transfer accuracy

Torque-safety feature prevents
damage to implant interface

Section off the groove
to shorten the abutment by 2mm

Two-Piece fixture-mount
with preparable abutment

Anatomical Profile ideal
for treatment in esthetic zone

Neck Matched to Major Diameter
Seals opening at crest of ridge
reducing need for bone grafting




www.implantdirect.eu | 00800 4030 4030



* Registered Trademarks of Nobel Biocare™

Features & Benefits

All-in-One Packaging InterActive™



- Ø: 3.2 3.7 4.3 5.0
- Lengths: 6 8 10 11.5 13 16
- Internal conical connection compatible with NobelActive™ and NobelReplace® CC*
- Revolutionary new fixture-mount that provides the accuracy of an open-tray transfer with the simplicity of a closed-tray transfer
- All-in-One Packaging for 160 Euro: implant, cover screw, healing collar, transfer and abutment
- Platform Switching to increase soft tissue volume
- Surface options: SBM textured or SBActive™ light HA coating

Drill Sequences

3.0mmD Platform 3.2mmD Implant			3.0mmD Platform 3.7mmD Implant			3.4mmD Platform 4.3mmD Implant			3.4mmD Platform 5.0mmD Implant		
		Optional 			Optional 			Optional 			Optional
SD2.3 2.3/2.0mmD	SD2.8 2.8/2.3mmD	CBD3.2	SD2.8 2.3/2.8mmD	SD3.4 3.4/2.8mmD	CBD3.7	SD3.4 3.4/2.8mmD or length specific NobelReplace® RP drill	SD3.8 3.8/3.4mmD or length specific NobelReplace® RP drill	CBD4.2	SD3.8 3.8/3.4mmD or length specific NobelReplace® RP drill	SD4.4 4.4/3.8mmD or length specific NobelReplace® RP drill	CBD4.7

Broaden your Prosthetic Options



Watch our online tutorial